

NO IMAGE
AVAILABLE

Pressure Injectable Arrowgard Blue Plus® Four-Lumen CVC

PI CVC KIT: 4-LUMEN 8.5 FR X 20 CM

SKU / Article #: EU-45854-HPSK

Highlights

Antimicrobial: Arrowgard Blue Plus*
Sharps Safety
CE Mark

Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Plus® Four-Lumen CVC

SKU / Article #: EU-45854-HPSK

Quantity	Component Description
1	Arrowg+ard Blue Plus® ¹ Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm, Pressure Injectable
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer with ECG Mark
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Injection Needle: Eclipse™ ² 25 Ga. x 1" (2.54 cm)
1	Pressure Transduction Probe
1	Syringe: 3 mL Luer-Lock
1	Syringe: 10 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 9 Fr. x 10.2 cm
4	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Drape: 47-1/4" x 59" (120 cm x 150 cm) with 4" (10 cm) fenestration, with adhesive
1	Safety Scalpel: #11
3	Prep Sponge Swab
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® ³ I.V. Advanced 3-1/4" x 4-1/2" (8.5 cm x 11.5 cm)
1	Dressing: GRIP-LOK® ⁴ Catheter Fixation Device


¹ Licensed under US Patent No. 7,329,412.

² A trademark of Becton, Dickinson and Company.

³ A registered trademark of 3M Company.

⁴ A registered trademark of TIDI.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.