

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Set

ARTERIAL CATH SET:18 GA CATH/20 GA NDL

SKU / Article #: FA-04018

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 25

Arterial Catheterization Set SKU / Article #: FA-04018	
Quantity	Component Description
1	Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 9-21/32" (24.5 cm)
1	Catheter: 18 Ga. x 4-1/4" (10.8 cm) Radiopaque over 20 Ga. TW Introducer Needle
1	Wing Clip: Suture

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.