

## Radial Artery Catheterization Set

### RA CATH SET: 20 GA X 3-1/16"

SKU / Article #: GH-04120-E

#### Highlights

CE Mark

#### Considerations


Not made with natural rubber latex.

**Sales Quantity/Case: 10**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Radial Artery Catheterization Set SKU / Article #: GH-04120-E	
Quantity	Component Description
1	Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
1	Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.