

Radial Artery Catheterization Set

RA CATH SET: 22 GA X 1-3/8 IN

SKU / Article #: GH-04122

Highlights

CE Mark

Considerations



Not made with natural rubber latex.

Sales Quantity/Case: 25

Representative photo(s): actual product specifications may vary from that shown in the photo

Radial Artery Catheterization Set SKU / Article #: GH-04122	
Quantity	Component Description
1	Catheter: 22 Ga. x 3.49 cm Radiopaque Polyurethane with Integral Suture Wing
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

