

Femoral Artery Catheterization Set

ARTERIAL CATH SET: 20 GA X 8IN

SKU / Article #: GH-04150

Highlights

CE Mark

Considerations



Not made with natural rubber latex.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Femoral Artery Catheterization Set SKU / Article #: GH-04150	
Quantity	Component Description
1	Catheter: 20 Ga. x 20 cm Radiopaque Polyurethane
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW

Documents



Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.