

HANDS-OFF® Two-Lumen CVC for use only with Arrow One-Piece PSI System

CVC COMPONENT: 2-LUMEN 7 FR X 30 CM

SKU / Article #: HO-14702

Highlights

CE Mark



Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

HANDS-OFF® Two-Lumen CVC for use only with Arrow One-Piece PSI System	
SKU / Article #: HO-14702	
Quantity	Component Description
1	Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm, Cath-Gard® Catheter Contamination Shield
2	Dust Cap: Non-Vented

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.