

NO IMAGE
AVAILABLE

Pressure Injectable Three-Lumen PICC
pre-loaded with Arrow® NaviCurve™
Stylet

PI PICC NAVICURVE: 3L 6FR x 55CM

SKU / Article #: IBAS-35563-CURV

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

†††MR Safe

Sales Quantity/Case: 5

Pressure Injectable Three-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet SKU / Article #: IBAS-35563-CURV	
Quantity	Component Description
1	Three Lumen TaperFree® Catheter†††: 6 Fr. (2.14 mm OD) x 55 cm, Pressure Injectable, Arrow® NaviCurve™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp


Pressure Injectable Three-Lumen PICC pre-loaded with Arrow NaviCurve® Stylet

SKU / Article #: IBAS-35563-CURV

Quantity	Component Description
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Sterile Procedure Sign
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
1	Transducer Cover with Gel and Rubber Bands
1	Remote Cover
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device
1	T-Piece Cover

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

