

NO IMAGE
AVAILABLE

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet

AGBA PICC NAVICURVE: 2L 5.5FR x 55CM

SKU / Article #: IBAS-45552-CURV

Highlights

Antimicrobial: Arrowg+ard Blue Advance*
Antithrombogenic: Arrowg+ard Blue Advance*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
‡‡‡MR Safe

Sales Quantity/Case: 5

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet SKU / Article #: IBAS-45552-CURV	
Quantity	Component Description
1	Two Lumen TaperFree® Catheter‡‡‡ with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 55 cm, Pressure Injectable, Arrow® NaviCurve™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock


Pressure Injectable Arrowgard Blue Advance Two-Lumen PICC pre-loaded with Arrow NaviCurve Stylet
SKU / Article #: IBAS-45552-CURV

Quantity	Component Description
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Sterile Procedure Sign
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
1	Transducer Cover with Gel and Rubber Bands
1	Remote Cover
1	Dressing: STATLOCK®2 PICC Plus Catheter Stabilization Device
1	T-Piece Cover

¹ Licensed under US Patent No. 7,329,412 and 11,219,706.

² A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.