



PSI Set for use with 7 - 7.5 Fr. Catheters

PSI SET: 8.5 FR

SKU / Article #: IC-09880

Considerations



Caution: This product contains dry natural rubber.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

PSI Set for use with 7 - 7.5 Fr. Catheters	
SKU / Article #: IC-09880	
Quantity	Component Description
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
1	Dust Cap: Non-Vented
1	Safety Scalpel: #11
1	Stopcock: 3-way High-Flow
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.