



## Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet

AGBA PICC NaviCurve: 1L 4.5FR x 55CM

SKU / Article #: IDLX-45541-CURV

### Highlights

- Antimicrobial: Arrowg+ard Blue Advance\*
- Antithrombogenic: Arrowg+ard Blue Advance\*
- Sharps Safety

### Considerations

- \*Chlorhexidine
- Not made with natural rubber latex.
- ##MR Safe

**Sales Quantity/Case: 3**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet

SKU / Article #: IDLX-45541-CURV

Quantity	Component Description
1	One Lumen TaperFree® Catheter## with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 4.5 Fr. (1.58 mm OD) x 55 cm, Pressure Injectable, Arrow® NaviCurve™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: Eclipse™ <sup>2</sup> 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow NaviCurve Stylet  
 SKU / Article #: IDLX-45541-CURV

Quantity	Component Description
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover with Gel and Rubber Bands
1	Remote Cover
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK® <sup>3</sup> PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow NaviCurve Stylet  
 SKU / Article #: IDLX-45541-CURV

Quantity	Component Description
1	Surgical Apparel: Mask
1	Tape: Steri-Strip®4
1	Valve: MicroCLAVE™5 Neutral Displacement Connector
1	T-Piece Cover

<sup>1</sup> Licensed under US Patent No. 7,329,412 and 11,219,706.


<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of C. R. Bard, Inc.

<sup>4</sup> A registered trademark of 3M Company.

<sup>5</sup> A trademark of ICU Medical, Inc.

## Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.