



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Two-Lumen Midline
PI MIDLINE 2L: 5FR X 20CM

SKU / Article #: IMID-32052-MAX

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
###MR Safe

Sales Quantity/Case: 5

Pressure Injectable Two-Lumen Midline SKU / Article #: IMID-32052-MAX	
Quantity	Component Description
1	Two-Lumen TaperFree® Catheter###: 5 Fr.(1.81 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with Extension Line Clamps and Contamination Guard
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener

Pressure Injectable Two-Lumen Midline

SKU / Article #: IMID-32052-MAX

Quantity	Component Description
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Overbed Sign
1	Patient ID Card
1	Checklist/CLIP Sheet
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover with Gel and Rubber Bands
1	Surgical Apparel: Bouffant Cap
1	Dressing: STATLOCK® ² PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® ³
2	Valve: MicroCLAVE™ ⁴ Neutral Displacement Connector

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of C. R. Bard, Inc.

³ A registered trademark of 3M Company.
⁴ A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.