



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen Midline

AGBA PI MIDLINE 2-L: 5.5FR X 15CM

SKU / Article #: IMID-41552-MAX

Highlights

Antimicrobial: Arrowg+ard Blue Advance*
Antithrombogenic: Arrowg+ard Blue Advance*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
†††MR Safe

Sales Quantity/Case: 5

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen Midline SKU / Article #: IMID-41552-MAX

Quantity	Component Description
1	Two-Lumen TaperFree® Catheter††† with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 15 cm Pressure Injectable with Blue FlexTip® Feature and Contamination Guard
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: Eclipse™ ² 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock

Pressure Injectable Arrowg+ard Blue Advance Two-Lumen Midline
SKU / Article #: IMID-41552-MAX

Quantity	Component Description
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Overbed Sign
1	Patient ID Card
1	Checklist/CLIP Sheet
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover with Gel and Rubber Bands
1	Surgical Apparel: Bouffant Cap
1	Dressing: STATLOCK® ³ PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® ⁴

Pressure Injectable Arrowg+ard Blue Advance Two-Lumen Midline
SKU / Article #: IMID-41552-MAX

Quantity	Component Description
2	Valve: MicroCLAVE™ ⁵ Neutral Displacement Connector

¹ Licensed under US Patent No. 7,329,412 and 11,219,706.
² A trademark of Becton, Dickinson and Company.
³ A registered trademark of C. R. Bard, Inc.
⁴ A registered trademark of 3M Company.
⁵ A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.