



Twin-Cath® 20/22 Two-Lumen Peripheral Catheter

2L INTRAVENOUS SET

SKU / Article #: IV-01100

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 25



Representative photo(s): actual product specifications may vary from that shown in the photo

Twin-Cath® 20/22 Two-Lumen Peripheral Catheter

SKU / Article #: IV-01100

Quantity	Component Description
1	Catheter: 18 Ga. x 1-3/4" (4.5 cm) Radiopaque over 22 Ga. RW Needle
1	Dust Cap: Non-Vented

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.