


Home / Central Venous Catheter / **JRT-43541-HPHNM**



Pressure Injectable Arrowg+ard Blue Advance® One-Lumen JACC with 80 cm Hydrophilic-Coated Nitinol Guidewire

JACC IR: 1-L 4.5FR NON-TIP 80 HMIT WIRE

SKU / Article #: JRT-43541-HPHNM

### Highlights

Antimicrobial: Arrowg+ard Blue Advance\*


Antithrombogenic: Arrowg+ard Blue Advance\*

### Considerations

\*Chlorhexidine

Not made with natural rubber latex.

**Sales Quantity/Case: 5**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Pressure Injectable Arrowg+ard Blue Advance® One-Lumen JACC with 80 cm Hydrophilic-Coated Nitinol Guidewire



SKU / Article #: JRT-43541-HPHNM

Quantity	Component Description
1	One-Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 4.5 Fr. (1.55 mm OD) x 35 cm, Pressure Injectable
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 4" (10 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	Catheter Trimmer
1	Safety Scalpel: #11
1	TightTrack™ Tuner: 9" Stainless Steel
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Paper Tape Measure
1	Dressing: Tegaderm® <sup>2</sup> 10 cm x 12 cm

<sup>1</sup> Licensed under US Patent Nos. 6,706,024, 6,872,195 and 7,329,412.

<sup>2</sup> A registered trademark of 3M Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.