

## Two-Lumen CVC for use only with Arrow MAC™ Two-Lumen Central Venous Access Device

**CVC COMPONENT: 2-LUMEN 7 FR**

SKU / Article #: ML-00702

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.



**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Two-Lumen CVC for use only with Arrow MAC™ Two-Lumen Central Venous Access Device  
SKU / Article #: ML-00702

Quantity	Component Description
1	Two-Lumen Catheter: 7 Fr. (2.5 mm OD) Catheter Contamination Guard
2	Dust Cap: Non-Vented

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.