



## One-Lumen CVC

CVC SET: S-L 16GA X 20CM

Microintroducer

SKU / Article #: MP-04301-CJ

### Considerations



Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

One-Lumen CVC	
SKU / Article #: MP-04301-CJ	
Quantity	Component Description
1	One-Lumen Catheter: 16 Ga. (1.8 mm OD) x 20 cm, Stylet
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 23-5/8" (60 cm) (Straight Stiff Tip on One End - "J" Tip on Other) with Arrow Advancer and hub
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 2-3/4" (7 cm) XTW
1	Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 1-3/8" (3.5 cm) XTW
1	Syringe: 2 mL Luer-Slip
1	Tissue Dilator: 6 Fr x 6.5 cm
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Y-hub
2	Guidewire Introducer

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.