



Two-Lumen CVC

CVC SET: 2-L 7FR X 30CM Microintroducer

SKU / Article #: MP-14702-CA

Contraindications & Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen CVC SKU / Article #: MP-14702-CA	
Quantity	Component Description
1	Two-Lumen Indwelling Catheter: 7 Fr x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps and Stylet
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 27-1/2" (70 cm) (Angled Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer and hub
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 2-3/4" (7 cm) XTW
1	Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 1-3/8" (3.5 cm) XTW
1	Syringe: 2 mL Luer-Slip
1	Tissue Dilator: 8 Fr x 9.7 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
3	Towel: 13" x 18" (33 cm x 46 cm)
1	Y-hub
2	Guidewire Introducer

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.