


Pressure Injectable Three-Lumen CVC PI CVC KIT: 3-L 7FR X 20CM Microintroduc

SKU / Article #: MPH-15703-JA

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Three-Lumen CVC SKU / Article #: MPH-15703-JA

Quantity	Component Description
1	Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Pressure Injectable, Stylet
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 23-5/8" (60 cm) (Angled Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer and hub
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 2-3/4" (7 cm) XTW
1	Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 1-3/8" (3.5 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Syringe: 10 mL Luer-Slip
1	Syringe: 2 mL Luer-Slip
1	Tissue Dilator: 8 Fr x 9.7 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Medicine Cup: 60 mL
1	SharpsAway® Disposal Cup



Pressure Injectable Three-Lumen CVC

SKU / Article #: MPH-15703-JA

Quantity	Component Description
1	Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive
1	Needle Holder
1	Scalpel: #11
2	Stopcock: 4-way
1	Prep Sponge Swab
6	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: OpSite ^{®1} 4-3/4" x 4" (12 cm x 10 cm)
1	Suture: 3-0 Silk with Curved Needle
1	Tray: Prep
1	Y-hub
2	Guidewire Introducer

¹ A registered trademark of T.J. Smith & Nephew, Limited.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.