

## Three-Lumen CVC

CVC KIT: 3-LUMEN 7FR X 16CM

SKU / Article #: MTO-12703-J

### Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5



*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Three-Lumen CVC

SKU / Article #: MTO-12703-J

Quantity	Component Description
1	Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Injection Needle: 23 Ga. x 1-1/4" (3.20 cm)
2	Syringe: 3 mL Luer-Slip
1	Syringe: 5 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8 Fr. (2.7 mm) x 9 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup
1	Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive
1	Needle Holder
1	Scalpel: #11
1	Dressing: Adhesive 5" x 4"(12.5 cm x 10 cm)
2	Prep Sponge Swab
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 2-0 Silk with Curved Needle
1	Tray: Prep

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.