

## Radial Artery Catheterization

### RA CATH KIT: 20 GA X 1-3/4"

SKU / Article #: NA-04020-K

#### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 10**






*Representative photo(s): actual product specifications may vary from that shown in the photo*

Radial Artery Catheterization SKU / Article #: NA-04020-K	
Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Syringe: 3 mL Luer-Lock
1	3 mL Applicator 2% CHG and 70% IPA ChloroPrep® <sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
1	Wing Clip: Suture
1	SharpsAway® II Locking Disposal Cup
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
1	Scalpel: #11
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.