

Representative photo(s): actual product specifications may vary from that shown in the photo

MAC™ Two-Lumen Central Venous Access for use with 7.5 - 8 Fr. Catheters
MULTI-LUMEN/PSI KIT

SKU / Article #: NA-11242-CDC

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
Includes Lidocaine
±Precaution: Contains Phthalate: DEHP


Sales Quantity/Case: 5

MAC™ Two-Lumen Central Venous Access for use with 7.5 - 8 Fr. Catheters SKU / Article #: NA-11242-CDC	
Quantity	Component Description
1	Two-Lumen Central Venous Access Device: 115 mm Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: SafetyGlide™1 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator

MAC Two-Lumen Central Venous Access for use with 7.5 - 8 Fr. Catheters SKU / Article #: NA-11242-CDC	
Quantity	Component Description
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
1	Dust Cap: Vented
2	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Needle Holder
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
1	Medication Label: 1% Lidocaine
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Silk with Curved Needle
1	Tubing±: 20" (50 cm) Extension

¹ A trademark of Becton, Dickinson and Company.
² A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.