



Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization

ARTERIAL CATH KIT: 20 GA X 5IN

SKU / Article #: NS-04100-S

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 5

| Arterial Catheterization SKU / Article #: NS-04100-S | |
|---|---|
| Quantity | Component Description |
| 1 | Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP |
| 1 | Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends) |
| 1 | Catheter: 20 Ga. x 1-1/4" (3.2 cm) Jelco® ¹ Catheter over 22 Ga. Introducer Needle |
| 1 | Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) |
| 1 | Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW |
| 2 | Syringe: 3 mL Luer-Lock |
| 1 | Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration |
| 1 | Prep Sponge Swab |
| 3 | Gauze Pad: 2" x 2" (5 cm x 5 cm) |

¹ A registered trademark of Medex, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.