



## Insertion Set for use with Arrow® 3 Fr. PICC

### ACCESS TRAY

SKU / Article #: PI-01351-LS

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 10**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Insertion Set for use with Arrow® 3 Fr. PICC

SKU / Article #: PI-01351-LS

Quantity	Component Description
1	GlideThru™ Peel-Away Sheath: 3 Fr. x 2-3/4" (7 cm) Radiopaque over 3 Fr. Dilator
1	Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On Both Ends)
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.