

## One-Lumen PICC with 80cm Wire

PICC SET: 3 FR X 21-1/2" (55 CM)

SKU / Article #: PR-01351-MH

### Highlights

CE Mark

### Considerations




Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

One-Lumen PICC with 80cm Wire SKU / Article #: PR-01351-MH	
Quantity	Component Description
1	One Lumen Catheter: 3 Fr. (1.03 mm OD) x 55 cm, Contamination Guard
1	GlideThru™ Peel-Away Sheath: 3 Fr. x 2-3/4" (7 cm) Radiopaque over 3 Fr. Dilator
1	Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On Both Ends)
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 31-1/2" (80 cm)(Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	Catheter Trimmer
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Paper Tape Measure

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.