

## Pressure Injectable One-Lumen PICC with 80 cm Wire

PICC SET: 4 FR X 19-5/8" (50 CM)

SKU / Article #: PR-05041-HP

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 5**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

### Pressure Injectable One-Lumen PICC with 80 cm Wire

SKU / Article #: PR-05041-HP

Quantity	Component Description
1	One Lumen Catheter: 4 Fr. (1.45 mm OD) x 50 cm, Pressure Injectable, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 31-1/2" (80 cm)(Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	PICC Line Insertion Checklist
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Paper Tape Measure
1	Dressing: STATLOCK®1 PICC Plus Catheter Stabilization Device

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.