

One-Lumen PICC with 45cm Wire

PICC SET: 4 FR X 19-1/2" (50 CM)

SKU / Article #: PR-05041

Highlights

CE Mark

Considerations




Not made with natural rubber latex.

Sales Quantity/Case: 5

One-Lumen PICC with 45cm Wire SKU / Article #: PR-05041	
Quantity	Component Description
1	One Lumen Catheter: 4 Fr. (1.45 mm OD) x 50 cm, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	Chart Sticker
1	Patient Information Booklet
1	Patient ID Card
1	Paper Tape Measure
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.