



Pressure Injectable Two-Lumen PICC with 80 cm Wire

PICC SET 2-LUMEN: 5 FR X 19-5/8" (50
CM)

SKU / Article #: PR-05052-HP

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

*Representative photo(s): actual product
specifications may vary from that shown in the
photo*



Pressure Injectable Two-Lumen PICC with 80 cm Wire

SKU / Article #: PR-05052-HP

Quantity	Component Description
1	Two Lumen Catheter: 5 Fr. (1.81 mm OD) x 50 cm, Pressure Injectable, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 4" (10 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 31-1/2" (80 cm)(Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	PICC Line Insertion Checklist
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Paper Tape Measure
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

