

## Two-Lumen PICC with 80cm Wire

PICC SET: 2-LUMEN 5 FR X 19-1/2" (50 CM)

SKU / Article #: PR-05052-MW

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5



*Representative photo(s): actual product specifications may vary from that shown in the photo*

Two-Lumen PICC with 80cm Wire  
SKU / Article #: PR-05052-MW

Quantity	Component Description
1	Two Lumen Catheter: 16 Ga. (1.76 mm OD) x 50 cm, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 4" (10 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 31-1/2" (80 cm)(Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Chart Sticker
1	Patient Information Booklet
1	Patient ID Card
1	Paper Tape Measure
1	Dressing: STATLOCK® <sup>1</sup> PICC Plus Catheter Stabilization Device

<sup>1</sup> A registered trademark of C. R. Bard, Inc.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

