

Two-Lumen PICC

PICC SET: 2-LUMEN 5 FR X 19-1/2" (50 CM)

SKU / Article #: PR-05052

Highlights

CE Mark

Considerations



Not made with natural rubber latex.

Sales Quantity/Case: 5

Two-Lumen PICC SKU / Article #: PR-05052	
Quantity	Component Description
1	Two Lumen Catheter: 16 Ga. (1.76 mm OD) x 50 cm, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 4" (10 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Chart Sticker
1	Patient Information Booklet
1	Patient ID Card
1	Paper Tape Measure
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.