

Pressure Injectable One-Lumen PICC with 80 cm Wire

PI PICC 1L 4FRX40CM W/80CM HYDRO
NIT SWG

SKU / Article #: PR-34041-HPHNM

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo


Pressure Injectable One-Lumen PICC with 80 cm Wire

SKU / Article #: PR-34041-HPHNM

Quantity	Component Description
1	One Lumen TaperFree [®] Catheter: 4 Fr. (1.45 mm OD) x 40 cm, Pressure Injectable, Blue FlexTip [®] and Contamination Guard
1	GlideThru [™] Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	SecondSite [™] Adjustable Hub: Fastener
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Checklist/CLIP Sheet
1	Paper Tape Measure
1	Dressing: STATLOCK ^{®1} PICC Plus Catheter Stabilization Device

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.