




Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Three-Lumen PICC with 80 cm Wire

PI PICC 3L 6FRX40CM W/80CM HYDRO NIT SWG

SKU / Article #: PR-34063-HPHNM

Highlights

CE Mark

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 5

Pressure Injectable Three-Lumen PICC with 80 cm Wire SKU / Article #: PR-34063-HPHNM	
Quantity	Component Description
1	Three Lumen TaperFree® Catheter: 6 Fr. (2.14 mm OD) x 40 cm, Pressure Injectable, Blue FlexTip® and Contamination Guard
1	GlideThru™ Peel-Away Sheath: 6 Fr. x 4" (10 cm) Radiopaque over 6 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Checklist/CLIP Sheet
1	Paper Tape Measure
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.