



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Two-Lumen PICC with 80 cm Wire

PI PICC 2L 5FRX55CM W/80CM HYDRO NIT SWG

SKU / Article #: PR-35552-HPHNM

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5






Pressure Injectable Two-Lumen PICC with 80 cm Wire

SKU / Article #: PR-35552-HPHNM

| Quantity | Component Description |
|----------|---|
| 1 | Two Lumen TaperFree® Catheter: 5 Fr. (1.81 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip® and Contamination Guard |
| 1 | GlideThru™ Peel-Away Sheath: 5 Fr. x 4" (10 cm) Radiopaque over 5 Fr. Dilator |
| 1 | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) |
| 1 | Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120) |
| 1 | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW |
| 1 | Syringe: 10 mL Luer-Lock |
| 2 | Dust Cap: Non-Vented |
| 1 | SecondSite™ Adjustable Hub: Fastener |
| 1 | SecondSite™ Adjustable Hub: Catheter Clamp |
| 1 | Catheter Trimmer |
| 1 | Safety Scalpel: #11 |
| 1 | Patient ID Card |
| 1 | Chart Sticker |
| 1 | Patient Information Booklet |
| 1 | Checklist/CLIP Sheet |
| 1 | Paper Tape Measure |
| 1 | Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device |

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.