

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen Midline

AGBA PI MIDLINE 2L: 5.5FR X 15CM

SKU / Article #: PR-41552-BAS

Highlights

Antimicrobial: Arrowg+ard Blue Advance*

Antithrombogenic: Arrowg+ard Blue Advance*

Sharps Safety

Considerations


*Chlorhexidine

Not made with natural rubber latex.

†††MR Safe

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo



Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen Midline


SKU / Article #: PR-41552-BAS

Quantity	Component Description
1	Two-Lumen TaperFree® Catheter††† with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 15 cm Pressure Injectable with Blue FlexTip® Feature and Contamination Guard
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Overbed Sign
1	Patient Information Booklet
1	Paper Tape Measure
1	Dressing: STATLOCK® ² PICC Plus Catheter Stabilization Device

¹ Licensed under US Patent No. 7,329,412 and 11,219,706.

² A registered trademark of C. R. Bard, Inc.

Documents



SCROLL TO VIEW MORE ➔

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.