

Pressure Injectable Arrowgard Blue Advance® Two-Lumen PICC

PI CG+ PICC SET: 2-L 5.5 FR X 40CM

SKU / Article #: PR-44052-BAS

Highlights

Antimicrobial: Arrowgard Blue Advance*
Antithrombogenic: Arrowgard Blue Advance*
Sharps Safety
CE Mark

Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC





SKU / Article #: PR-44052-BAS

Quantity	Component Description
1	Two Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 40 cm, Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Dressing: STATLOCK® ² Catheter Stabilization Device

¹ Licensed under US Patent No. 7,329,412.

² A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.