




Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC

PI AGBA PICC SET: 1-L 4.5 FR X 45CM

SKU / Article #: PR-44541-BAS

Highlights


Antimicrobial: Arrowg+ard Blue Advance*
 Antithrombogenic: Arrowg+ard Blue Advance*
 Sharps Safety
 CE Mark

Considerations

*Chlorhexidine
 Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo




Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC SKU / Article #: PR-44541-BAS	
Quantity	Component Description
1	One Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 4.5 Fr. (1.58 mm OD) x 45 cm, Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Dressing: STATLOCK® ² Catheter Stabilization Device

¹ Licensed under US Patent No. 7,329,412.

² A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.