




Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC

PI CG+ PICC SET: 2-L 5.5 FR X 55CM

SKU / Article #: PR-45552-BAS

Highlights

Antimicrobial: Arrowg+ard Blue Advance*
 Antithrombogenic: Arrowg+ard Blue Advance*
 Sharps Safety
 CE Mark

Considerations

*Chlorhexidine
 Not made with natural rubber latex.

Sales Quantity/Case: 5

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC


SKU / Article #: PR-45552-BAS

Quantity	Component Description
1	Two Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 55 cm, Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Dressing: STATLOCK® ² PICC Plus Catheter Stabilization Device

¹ Licensed under US Patent No. 7,329,412.

² A registered trademark of C. R. Bard, Inc.

Documents



SCROLL TO VIEW MORE ➡

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.