

## Pressure Injectable Arrowgard Blue Advance® Two-Lumen PICC

AM PICC SET 2-L: 5.5 FR X 55 CM

SKU / Article #: PR-45552-HP

### Highlights

---

Antimicrobial: Arrowgard Blue Advance\*

Antithrombogenic: Arrowgard Blue Advance\*

### Considerations

---

\*Chlorhexidine

Not made with natural rubber latex.

**Sales Quantity/Case: 5**

---

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Pressure Injectable Arrowgard Blue Advance® Two-Lumen PICC


SKU / Article #: PR-45552-HP

Quantity	Component Description
1	Two Lumen TaperFree® Catheter with Arrowgard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 5.5 Fr. (1.88 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip®
1	Peel-Away Sheath: 6 Fr. x 4" (10 cm) Radiopaque over 6 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Checklist/CLIP Sheet
1	Paper Tape Measure
1	Dressing: STATLOCK® <sup>2</sup> PICC Plus Catheter Stabilization Device

<sup>1</sup> Licensed under US Patent No. 7,329,412.

<sup>2</sup> A registered trademark of C. R. Bard, Inc.

## Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

