

One-Lumen PICC

PICC SET: 1-L 4 FR X 19-1/2" (50 CM)

SKU / Article #: PS-01451-NM

Considerations

Not made with natural rubber latex.


Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

One-Lumen PICC SKU / Article #: PS-01451-NM	
Quantity	Component Description
1	One Lumen Catheter: 4 Fr. (1.45 mm OD) x 50 cm, Blue FlexTip®
1	Peel-Away Sheath: 15 Ga. x 1-1/4" (3.18 cm) Radiopaque over 17 Ga. TW Introducer Needle
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.