

## One-Lumen PICC

PICC SET: 1L 16 GA X 27 1/2IN (70 CM)

SKU / Article #: PS-01671

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.




**Sales Quantity/Case: 10**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

One-Lumen PICC SKU / Article #: PS-01671	
Quantity	Component Description
1	One Lumen Catheter: 16 Ga. (1.80 mm OD) x 70 cm, Side Holes, Blue FlexTip®, Contamination Guard
1	Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Patient Information Booklet
1	Patient ID Card
1	Dressing: STATLOCK® <sup>1</sup> PICC Plus Catheter Stabilization Device

<sup>1</sup> A registered trademark of C. R. Bard, Inc.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.