

Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen PICC
PICC SET: 2-L 16 GA X 27-1/2 IN (70 CM)

SKU / Article #: PS-01672

Highlights

CE Mark

Considerations

Not made with natural rubber latex.


Sales Quantity/Case: 10

Two-Lumen PICC
SKU / Article #: PS-01672

Quantity	Component Description
1	Two Lumen Catheter: 16 Ga. (1.76 mm OD) x 70 cm, Blue FlexTip®, Contamination Guard
1	Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
2	Dust Cap: Non-Vented
2	Clamp: Extension Line
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Patient Information Booklet
1	Patient ID Card
1	Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.