

NO IMAGE  
AVAILABLE

Spring-Wire Guide .018" (.46 mm)  
Diameter x 31-1/2" (80 cm) Wire Proximal  
End: Solid Wire .018" Diameter Distal  
End: Coil .018" Diameter x 7 cm  
SPRING WIRE GUIDE COMPONENT

SKU / Article #: PW-18080

Highlights

CE Mark

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 5

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm  
SKU / Article #: PW-18080

| Quantity | Component Description  |
|----------|--|
| 1        | Spring-Wire Guide: .018" (0.46 mm) dia. x 31-1/2" (80 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.