

NO IMAGE
AVAILABLE

Spring-Wire Guide .018" (.46 mm)
Diameter x 51-3/16" (130 cm) Wire.
Proximal End: Solid Wire .018" Diameter
Distal End: Coil .018" Diameter x 7 cm
SPRING WIRE GUIDE COMPONENT

SKU / Article #: PW-18130

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Spring-Wire Guide .018" (.46 mm) Diameter x 51-3/16" (130 cm) Wire. Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm SKU / Article #: PW-18130	
Quantity	Component Description
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

