



Spring-Wire Guide .018" (.46 mm) Diameter x 51-3/16" (130 cm) Wire.  
 Proximal End: Solid Wire .018" Diameter  
 Distal End: Coil .018" Diameter x 7 cm  
**SPRING WIRE GUIDE COMPONENT**

SKU / Article #: PW-18130

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Spring-Wire Guide .018" (.46 mm) Diameter x 51-3/16" (130 cm) Wire. Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm  
 SKU / Article #: PW-18130

Quantity	Component Description
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

**Documents**



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

