



## Radial Artery Catheterization Set

RA CATH SET: 18 GA X 1 3/4IN

SKU / Article #: RA-04018

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 50**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Radial Artery Catheterization Set SKU / Article #: RA-04018	
Quantity	Component Description
1	Catheter: 18 Ga. x 4.45 cm Radiopaque FEP over 20 Ga. TW Introducer Needle with Integral .025" (0.64 mm) dia. Spring-Wire Guide
1	Wing Clip: Suture

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

