

Radial Artery Catheterization Set

RA CATH SET: 20 GA X 1-3/4IN

SKU / Article #: RA-04020

Highlights

CE Mark

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 50

Representative photo(s): actual product specifications may vary from that shown in the photo

Radial Artery Catheterization Set SKU / Article #: RA-04020	
Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Wing Clip: Suture

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

