



## Radial Artery Catheterization Set

RA CATH SET: 22 GA X 1 3/4IN (4.45 CM)

SKU / Article #: RA-04022

### Highlights

CE Mark

### Considerations


Not made with natural rubber latex.

**Sales Quantity/Case: 50**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Radial Artery Catheterization Set SKU / Article #: RA-04022	
Quantity	Component Description
1	Spring-Wire Guide/Tube Assembly: .018" (0.46 mm) dia. x 5-3/8" (14 cm)
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Wing Clip: Suture

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.