



Representative photo(s): actual product specifications may vary from that shown in the photo

Radial Artery Catheterization Set

RA CATH SET: 22 GA X 1 3/4IN (4.45 CM)

SKU / Article #: RA-04022

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 50

Radial Artery Catheterization Set SKU / Article #: RA-04022	
Quantity	Component Description
1	Spring-Wire Guide/Tube Assembly: .018" (0.46 mm) dia. x 5-3/8" (14 cm)
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Wing Clip: Suture

Documents



Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.