



## QuickFlash® Radial Artery Catheterization Set

RA CATH SET: 20 GA X 1 1/2IN

SKU / Article #: RA-04220-W

### Highlights

CE Mark

### Considerations


Not made with natural rubber latex.

**Sales Quantity/Case: 50**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

QuickFlash® Radial Artery Catheterization Set	
SKU / Article #: RA-04220-W	
Quantity	Component Description
1	Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.