



Representative photo(s): actual product specifications may vary from that shown in the photo

QuickFlash® Radial Artery Catheterization Set

RA CATH SET: 20 GA X 1-1/2 IN

SKU / Article #: RA-04220

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 50

| QuickFlash® Radial Artery Catheterization Set SKU / Article #: RA-04220 | |
|--|--|
| Quantity | Component Description |
| 1 | QuickFlash® Radial Artery Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.