



RIC® Rapid Infusion Catheter Exchange Set

TRAUMA SET: 8.5 FR

SKU / Article #: RC-09850

Highlights

CE Mark

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 25

Representative photo(s): actual product specifications may vary from that shown in the photo

| RIC® Rapid Infusion Catheter Exchange Set SKU / Article #: RC-09850 | |
|--|--|
| Quantity | Component Description |
| 1 | Sheath: 8.5 Fr. x 2-1/2" (6.4 cm) Arrow-Flex® Radiopaque FEP with Tissue Dilator |
| 1 | Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends) |
| 1 | Scalpel: #11 |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.